

EXPERIMENTAL PROTOCOL: Maintenance, passaging and cryopreservation of colon cancer organoids

Author/s: Mantrana A., Guil-Luna S., Aranda-Aguilar E., Rodríguez Ariza A.

Contact information: ana.mantrana@imibic.org; v22gulus@uco.es

Introduction

This protocol is for the maintenance, passaging and cryopreservation of colorectal cancer organoids.

Reagents

- Colorectal Cancer Organoid Media (see below)
- Dimethyl sulfoxide (DMSO) (Sigma-Aldrich, cat. no. 34943)
- Falcon tubes, 15 mL (Sarstedt, cat. no. 62.554.502)
- FBS (Biowest, S1810)
- Matrigel (Corning #354230)
- Organoid Culture Media (see below)
- TrypLE Express (Gibco, cat. no. 12604-013)

Instruments and tools

- Mr Frosty freezing container (ThermoFisher, 5100-0036)
- Nunc Cryovials (ThermoFisher, cat. no. 375418)
- Plates, 24 well (Greiner, cat. no. 662160)

Before start

Preparation of Reagents:

Component	Company/Cat#	Stock Conc.	Volume	Final Conc.
Advanced DMEM/F12	Sigma #D6421	500 mL	486.9 mL	1x
HEPES	Media Kitchen (Life Technologies #15630080)	1 M	5 mL	10 mM
Glutamax	Media Kitchen (Life Technologies #35050-061)	100x (200 mM)	5 mL	1x (2 mM)
Zell Shield	Minerva Biolabs GmbH, Berlin, Germany	100x	5 mL	1x

Colorectal Cancer Organoid Media (CRC Medium)

Component	Company/Cat#	Final Conc.
OC Medium	See above	1x
A 83-01	TOCRIS #2939-10mg	0.5 μ M
B27	ThermoFisher Sc. #12587010	2x
EGF	Sigma Aldrich #SRP3027-500UG	0.25 μ g/ml
Gastrin	Sigma Aldrich #G9145-0.5MG	1 μ g/mL (480.7461nM)
N acetyl Cyst	Sigma Aldrich #A9165-5G	1.23mM
SB202190	Sigma Aldrich #S7067-5MG	5.0296 μ M
Y27632	Sigma Aldrich #Y0503-5mg	10 μ M

PROCEDURE

Priorities:

- 1) Maintenance organoids
- 2) Expansion organoids
- 3) Production of cryopreserved organoids biobank

Maintenancance of Organoids (Media change every Monday and Thursday)

1. Carefully aspirate 500 μ L medium from each well containing organoids.
2. Add 500 μ L of pre-warmed (37°C) CRC medium to each well.

Passaging of organoids

Passaging organoids without dissociation:

1. Carefully remove culture medium from each well and add 400 μ L of cold PBS per well.
2. Scrap matrigel with p1000 and transfer organoids to a 15 mL falcon tube.
3. Wash each well again with 400 μ L of cold PBS and transfer remaining organoids to the 15 mL falcon tube.
4. With a 21g syringe, separate matrigel from organoids mechanically.
5. Centrifuge at 300g for 5 minutes at 4°C and remove supernatant.
6. Resuspend pelleted organoids in the desired matrigel volume (dilution from 1 well \Rightarrow 2 wells or from 1 well \Rightarrow 4 wells depending on the initial density obtained)
7. Allow matrigel to settle for 30 minutes at 37°C.
8. Add 500 μ L of colorectal cancer organoid medium per well.

Passaging organoids with Dissociation:

1. Carefully remove culture medium from each well and add 500 μ L of Tryple per well.
2. Scrap matrigel with p1000 and incubate for 15 minutes at 37°C.
3. Check dissociation every 5 minutes and pipette with p1000 to enhance cell detachment.
4. Collect dissociated organoids into a 15 mL falcon tube and top up to 15 mL with PBS.

5. Centrifuge at 300g for 5 minutes RT, and remove supernatant.
6. Resuspend pelleted organoids in the desired matrigel volume (dilution from 1 well \Rightarrow 2 wells or from 1 well \Rightarrow 4 wells depending on the initial density obtained)
7. Allow matrigel to settle for 30 minutes at 37°C.
8. Add 500 μ L of colorectal cancer organoid medium per well.

Cryopreservation of organoids

1. Starting from passage 2, set aside 2 wells from each passage for cryopreservation.
2. Scrap the matrigel at the bottom of each well and transfer to 15 mL Falcon tubes.
3. Rewash each well with 500 μ L of organoid culture medium.
4. Top remaining Falcon tube with organoid culture medium up to 10 mL.
5. Centrifuge at 1500 rpm for 5 minutes at 4°C
6. Aspirate supernatant to the top of matrigel
7. Add FBS with 10% DMSO (900 μ L FBS + 100 μ L DMSO), using 1 mL of this mix per 2 wells.
8. Gently resuspend organoids and transfer 1mL of freezer mix to each cryopreservation vial.
9. Place the vials into a Mr. Frosty and store at -80°C for at least 24h. Within several days, transfer vials to liquid nitrogen storage.